

# Application Sheet

## H. pylori IgG

**REF 51222**

## Instrument: Elisys Quattro

**REF 16300**

### Setting:

**HPYLG\_QT\_51222\_R4**

**HPYLG\_CO\_51222\_R4**

For all essential product information, please refer to the User Manual of the analyser and the Instructions for Use for reagents, controls and calibrators. Please refer to the Instructions for Use for performance data.

This Application Sheet provides additional information regarding the use of the assay on the HUMAN ELISA analyser.

The parameters defined in this application sheet have been developed to provide optimal product performance with the assay and instrument combination. Any modification to these parameters may negatively affect the performance and results of the assay. The user is responsible for the validation of any modification to the protocol here described.

### Material Required

Material	REF	Comments
H. pylori IgG	51222	
Reagent Rack 0 Elisys Quattro	16300/38	All reagents, except DIL
Sample Rack 20 Elisys Quattro	16300/41	Sample collection tubes
Diluent prepared in	17350/B100	Reagent bottle 100 ml

### Additional Notes

Due to technical reasons, reagent DIL must be transferred to a disposable bottle that fits into the reagent rack (Reagent Bottle 100 ml, REF 17350/B100).

Chose a setting according to the preferred option for data analysis: QT for quantitative results (U/ml) or CO for results in relation to the cut-off (Abs/COV).

OD-values of 3.5 or above are outside the photometer's dynamic range. For this reason, they are displayed as 4.0.