

Application Sheet

CMV IgM

REF 51103

Instrument: Elisis Quattro

REF 16300

Setting:

CMVM_51103_R4.asy

For all essential product information, please refer to the User Manual of the analyzer and the instructions for use for reagents, controls, calibrators and the Certificate of Analysis. Please refer to the instruction for use for performance data.

This Application Sheet provides additional information regarding the use of the assay on the HUMAN ELISA analyzer.

The parameters defined in this application sheet have been developed to provide optimal product performance with the assay and instrument combination. Any modification to these parameters may negatively affect the performance and results of the assay. The user is responsible for the validation of any modification to the protocol here described.

Material Required

Material	REF	Comments
CMV IgM	51103	
Reagent Rack 0 Elisis Quattro	16300/38	NC, PC, CON; SUB; STOP
Sample Rack 20 Elisis Quattro	16300/41	Suitable for sample collection tubes

Additional Notes

The second incubation time is displayed in the report as the sum of incubating the sample pre-dilution in the tube (target: 5 min) and the incubation in the plate (target: 30 min). Since the resulting time is compared with a specification of 30 min (plate incubation), the reported time is flagged in the report. The time log in the worklist allows to check that incubation took the correct time.