

Application Sheet

Rub IgG

REF 51208

Instrument: HumaReader Single plus

REF 18000

Setting:

RUBG_CO_51208_R2

RUBG_QT_51208_R6

For all essential product information, please refer to the User Manual of the analyzer and the instructions for use for reagents, controls, calibrators and the Certificate of Analysis. Please refer to the instruction for use for performance data.

This Application Sheet provides additional information regarding the use of the assay on the HUMAN ELISA analyzer.

The parameters defined in this application sheet have been developed to provide optimal product performance with the assay and instrument combination. Any modification to these parameters may negatively affect the performance and results of the assay. The user is responsible for the validation of any modification to the protocol here described.

Additional Notes

The results should be interpreted with regard to a cut-off value. Therefore use the cut-off assay file first.

The QT assay is programmed to get concentrations for positive samples in the report. For this reason all results are shown with concentrations. Concentrations below 10 U/ml and above 65 U/ml are calculated by extrapolation and should be considered as "out of range".

Regard that the reader needs special positioning of controls. Positioning differs for cut-off and for quantitative analysis. Note: 8-well-strips are breakable.

Blank – A1
CC – B1/C1
PC – D1/E1
NC – F1/G1

The validity criteria can't be checked by the reader. They have to be checked by the user.

Exception: Blank < 0,15.