

# Application Sheet

## RUBELLA IgG

**REF 51208**

## Instrument: Elisys Quattro

**REF 16300**

### Setting:

**RUBG\_51208\_R6**

For all essential product information, please refer to the User Manual of the analyzer and the instructions for use for reagents, controls, calibrators and the Certificate of Analysis. Please refer to the instruction for use for performance data.

This Application Sheet provides additional information regarding the use of the assay on the HUMAN ELISA analyzer.

The parameters defined in this application sheet have been developed to provide optimal product performance with the assay and instrument combination. Any modification to these parameters may negatively affect the performance and results of the assay. The user is responsible for the validation of any modification to the protocol here described.

### Material Required

Material	REF	Comments
Rubella IgG	51208	
Reagent Rack 0 Elisys Quattro	16300/38	NC, CC, PC, CON, SUB, STOP
Sample Rack 20 Elisys Quattro	16300/41	Suitable for sample collection tubes, Dilution tubes

### Additional Notes

Although the results should be interpreted with regard to a cut-off value, the assay is programmed as a quantitative assay to get concentrations for positive samples in the report. For this reason even some equivocal results are shown with concentrations.

Concentrations below 10 U/ml and above 65 U/ml are calculated by extrapolation and should be considered as “out of range”.