

# HUMAN

Target Values

## BioRad Lyphocheck Immunoassay Plus Control

370 Trilevel 12x5ml

Verison 03/ 11 - 2021

371 Level 1 12x5ml

Level 1 40401

**REF**


372 Level 2 12x5ml

**LOT**

40400

Level 2 40402

**IVD**

 2023-10-31

373 Level 3 12x5ml

Level 3 40403

370X MiniPak 3x5ml

Link for latest version of target value sheet:

<https://www.human.de/cli/target-value-sheets/>



<b>REF</b>	Assay Parameter Ensayo	Valid from Gültig ab Válido a partir de	Units Einheiten Unidades	Target Sollwert Valor asignado Cible	Range Bereich Intervalo	Target Sollwert Valor asignado Cible	Range Bereich Intervalo	Target Sollwert Valor asignado Cible	Range Bereich Intervalo
	Essai	Valable à partir de <b>RG</b> Lot	Unité	Intervalle	Intervalle	Intervalle	Intervalle	Intervalle	Intervalle
					Level 1 - 40401	Level 2 - 40402		Level 3 - 40403	
84600	TSH	20001	µIU/ml	0.68	0.45 - 0.92	6.08	4.00 - 8.25	35.14	23.10 - 47.69
84620	fT3	20001	pg/ml	2.49	1.64 - 3.38	6.46	4.25 - 8.77	9.59	6.30 - 13.02
84615	T3	20001	ng/ml	0.80	0.53 - 1.09	2.14	1.41 - 2.90	3.35	2.20 - 4.55
84610	fT4	20001	pg/ml	6.00	3.94 - 8.14	14.10	9.27 - 19.14	26.59	17.48 - 36.09
84605	T4	20001	ng/ml	53.66	35.27 - 72.82	112.6	73.99 - 152.8	154.29	101.42 - 209.39
82120	FSH	20001	mIU/ml	5.42	3.56 - 7.36	23.43	15.40 - 31.80	43.29	28.46 - 58.75
82115	LH	20001	mIU/ml	2.13	1.40 - 2.89	32.18	21.15 - 43.7	118.25	77.73 - 160.48
83130	PSA	20001	ng/ml	1.08	0.710 - 1.47	4.20	2.76 - 5.70	27.17	17.86 - 36.87
83125	fPSA	20001	ng/ml	0.36	0.24 - 0.49	1.96	1.29 - 2.66	11.85	7.79 - 16.08
83100	AFP	20001	ng/ml	16.61	9.97 - 23.25	87.8	52.68 - 122.92	191.93	115.16 - 268.70
82100	HCG	20001	mIU/ml	5.13	3.37 - 6.96	23.16	15.22 - 31.43	161.89	106.41 - 219.71
82105	Progesterone*	20001	ng/ml	<0,1	<0,1 - <0,1	12.69	8.34 - 17.22	33.04	21.72 - 44.84
82110	Testosterone*	20001	ng/ml	0.90	0.54 - 1.26	3.41	2.05 - 4.77	6.09	4.14 - 9.66

\*Conformity assessment pending