


HUMAN

Target Values

BioRad Lyphocheck Tumor Marker Plus Control

Version 01/ 05 - 2022

REF 367 Level 1 6x2ml Level 1 74601
LOT 368 Level 2 6x2ml 74600 Level 2 74602 IVD  2024-05-28
 369 Level 3 6x2ml Level 3 74603
 368X Trilevel MiniPak 3x2ml

Link for latest version of target value sheet:

<https://www.human.de/cia/target-value-sheets/>



REF	Assay Parameter Ensayo	Valid from Gültig ab Válido a partir de	Units Einheiten Unidades	Target Sollwert Valor asignado Cible	Range Bereich Intervalo	Target Sollwert Valor asignado Cible	Range Bereich Intervalo	Target Sollwert Valor asignado Cible	Range Bereich Intervalo
	Essai	Valable à partir de [RGT] Lot	Unité	Intervalle	Intervalle	Intervalle	Intervalle	Intervalle	
				Level 1 - 74601		Level 2 - 74602		Level 3 - 74603	
83130	PSA	20001	ng/ml	0.17	0.11 - 0.23	5.55	4.00 - 7.10	28.3	20.4 - 36.2
83125	fPSA	20001	ng/ml	0.10	0.07 - 0.13	3.07	2.30 - 3.84	14.4	10.8 - 18.0
83100	AFP	20001	ng/ml	10.3	7.84 - 12.8	80.2	60.9 - 99.4	228	173 - 283
83115	CA 15-3	21001	U/ml	22.3	15.6 - 29.0	59.4	41.6 - 77.2	124	86.7 - 161
83110	CA 125	21001	U/ml	36.9	27.7 - 46.1	95.8	71.9 - 120	283	212 - 353
83120	CA 19-9	21001	U/ml	14.7	10.4 - 18.9	31.6	22.5 - 40.8	90.8	64.5 - 117
83105	CEA	22001	ng/ml	2.56	1.92 - 3.20	19.0	14.2 - 23.7	63.1	47.3 - 78.9