

HUMAN

Target Values

BioRad Lyphocheck Immunoassay Plus Control

370 Trilevel 12x5ml
371 Level 1 12x5ml
372 Level 2 12x5ml
373 Level 3 12x5ml
370X MiniPak 3x5ml

Version 01/ 09-2023

REF

LOT 40430

Level 1 40431
Level 2 40432
Level 3 40433

IVD

2025-10-31

Link for latest version of target value sheet:

<https://www.human.de/cia/target-value-sheets/>



REF	Assay Parameter Ensayo Essai	Valid from Gültig ab Válido a partir de Valable à partir de [RGT] Lot	Units Einheiten Unidades Unité	Target		Range		Target		Range	
				Sollwert Valor asignado Cible	Valor asignado Cible	Bereich Intervalo Intervalle	Intervalo	Sollwert Valor asignado Cible	Valor asignado Cible	Bereich Intervalo Intervalle	Intervalle
				Level 1 - 40431		Level 2 - 40432		Level 3 - 40433			
84600	TSH	23001	µIU/ml	0.530	0.340 - 0.720	7.40	4.81 - 9.99	43.3	28.2 - 58.5		
84620	fT3	23002	pg/ml	2.90	1.89 - 3.92	7.56	4.91 - 10.2	10.5	6.82 - 14.2		
84615	T3	23001	ng/ml	0.900	0.590 - 1.22	2.43	1.58 - 3.28	3.54	2.30 - 4.78		
84610	fT4	23001	pg/ml	6.43	4.18 - 8.68	14.5	9.45 - 19.6	25.6	16.6 - 34.5		
84605	T4	23001	ng/ml	60.1	39.1 - 81.1	117	76.3 - 158	148	96.4 - 200		
82120	FSH	23001	mIU/ml	7.74	5.03 - 10.5	23.8	15.5 - 32.2	52.6	34.2 - 70.9		
82115	LH	23001	mIU/ml	1.84	1.20 - 2.48	23.9	15.5 - 32.2	91.1	59.2 - 123		
83130	PSA	23001	ng/ml	0.780	0.510 - 1.05	3.98	2.59 - 5.37	30.2	19.6 - 40.7		
83125	fPSA	23001	ng/ml	0.490	0.320 - 0.660	2.01	1.31 - 2.71	11.4	7.38 - 15.3		
83100	AFP	23001	ng/ml	19.1	12.4 - 25.8	108	70.1 - 146	233	151 - 314		
82100	HCG	23001	mIU/ml	5.16	3.35 - 6.97	21.7	14.1 - 29.4	209	136 - 283		
82105	Progesterone	23001	ng/ml	0.140	<0,100 - 0.190	12.7	8.26 - 17.2	28.3	18.4 - 38.3		
82110	Testosterone	23001	ng/ml	0.830	0.540 - 1.12	3.68	2.39 - 4.97	7.86	5.11 - 10.6		
85000	25-OH Vitamin D	23001	ng/ml	22.2	14.4 - 29.9	25.4	16.5 - 34.3	26.1	17.0 - 35.2		
83105	CEA	23001	ng/ml	2.80	1.82 - 3.78	18.7	12.2 - 25.3	35.2	22.9 - 47.5		
82130	Estradiol	23001	pg/ml	52.4	34.1 - 70.7	273	178 - 369	688	447 - 928		
85049	Ferritin*	23001	ng/ml	74.3	48.3 - 100	201	131 - 272	455	296 - 614		
82125	Prolactin	23001	ng/ml	9.06	5.89 - 12.2	24.3	15.8 - 32.9	50.1	32.6 - 67.6		
85055	Vitamin B12	23001	pg/ml	427	278 - 577	702	456 - 948	1045	679 - 1411		

* No CE - only for use outside Europe