

HC-CONTROL

Intended for / Verwendung für / para el uso con / pour l'usage avec
HUMACOUNT, HUMACOUNT Plus

REF 17400/40

Target Values / Sollwerte / Valores Meta / Valeurs Souhaitées

		LOT 192					
Parameter		LEVEL 1		LEVEL 2		LEVEL 3	
		low		normal		high	
		MW	± R	MW	± R	MW	± R
WBC	x 10 ⁹ cells/l	1.8	± 0.4	7.4	± 0.6	20.2	± 2.0
LYM	x 10 ⁹ cells/l	1.0	± 0.3	2.2	± 0.4	2.6	± 0.8
MID	x 10 ⁹ cells/l	0.2	± 0.1	0.6	± 0.2	0.9	± 0.4
GRA	x 10 ⁹ cells/l	0.6	± 0.2	4.6	± 0.6	16.6	± 1.2
LYM	%	55.4	± 5.0	29.7	± 4.0	13.0	± 4.0
MID	%	13.4	± 4.0	7.6	± 4.0	4.6	± 4.0
GRA	%	31.2	± 4.0	62.7	± 5.0	82.5	± 5.0
RBC	x 10 ¹² cells/l	2.57	± 0.12	4.87	± 0.15	5.74	± 0.20
HGB	g/dl	6.0	± 0.4	13.3	± 0.5	17.8	± 0.6
HGB	mmol/l	3.7	± 0.3	8.3	± 0.4	11.0	± 0.4
HCT	%	20.8	± 1.5	43.5	± 2.0	56.2	± 2.5
MCV	fl	81	± 4	89	± 4	98	± 4
MCH	pg	23.3	± 1.5	27.3	± 1.5	31.0	± 2.0
MCHC	g/dl	28.8	± 2.0	30.5	± 2.0	31.7	± 2.5
RDW-CV	%	17.9	± 3.0	16.8	± 3.0	15.8	± 3.0
PLT	x 10 ⁹ cells/l	62	± 15	272	± 30	540	± 55
PCT	%	0.05	± 0.02	0.22	± 0.03	0.43	± 0.04
MPV	fl	7.7	± 1.5	8.0	± 1.5	8.0	± 1.5
PDW-CV	%	39.0	± 2.5	37.6	± 2.5	35.8	± 3.0

MW: Mean value / Mittelwert / Valor media / Valeur moyen

R: Expected range / Erwarteter Bereich / Rango esperado / Domaine espéré

CS-HC3 INF 17400401 04-2010-31

Human

Human Gesellschaft für Biochemica und Diagnostica mbH
Max-Planck-Ring 21 - D-65205 Wiesbaden - Germany

Tel: +49 6122 9988 0 - Fax: +49 6122 9988 100 - eMail: human@human.de