



Accessories for LED Microscopes

REF	Description	Fluo	Premium	Advanced	Classic	Light
14000/72	Blue filter	X	X	X	X	X
14000/54	C-mount adapter 0.4x for trinocular head	X	X	X	X	
14000/55	C-mount adapter 1x for trinocular head		X	X	X	
14000/80	C-mount adapter 0.5x for trinocular	X	X	X	X	
16160	Cleaner set (lens cleaning fluid, lens cleaning rods, lens cleaning tissues, air duster microfiber cloth, dust bellow)	X	X	X	X	X
14000/52	Dark field condenser		X	X	X	X
14000/53	Dark field condenser oil	X	X	X	X	X
14000/16	Eyepiece 10x 18 mm view field					X
14000/17	Eyepiece 10x 18 mm view field			X	X	
14000/18	Eyepiece 10x 20 mm view field		X			
14000/19	Eyepiece 16x 11 mm view field		X	X	X	X
14000/20	Eyepiece 20x 11 mm view field		X	X	X	
14000/133	Eyepiece 20x 22 mm view field	X				
14000/205	Eyepiece adapter 0.5x	X	X	X	X	X
14000/125	Eyepiece pointer achromatic 10x 18 mm					X
14000/128	Eyepiece pointer achromatic 16x 11 mm			X	X	X
14000/131	Eyepiece pointer fluo 10x 22 mm	X				
14000/132	Eyepiece pointer fluo 15x 22 mm	X				
14000/127	Eyepiece pointer infinity 10x 20 mm		X			
14000/130	Eyepiece pointer infinity 10x 22 mm	X				
14000/21	Eyepiece micrometer 10 mm with 100 division for 18 mm view field			X	X	X
14000/22	Eyepiece micrometer 10 mm with 100 division for 20 mm view field		X			
14000/126	Eyepiece pointer plan achromatic 10x 18 mm			X	X	
14000/73	Frosted filter	X	X	X	X	X
14000/74	Green filter	X	X	X	X	X
14000/71	Mirror set			X	X	X
14000/10	Monocular head			X	X	
14000/23	Objective achromatic 4x					X
14000/24	Objective achromatic 10x					X
14000/30	Objective achromatic 20x					X
14000/25	Objective achromatic 40x					X
14000/26	Objective achromatic 60x					X
14000/27	Objective achromatic 100x					X
14000/32	Objective achromatic 60x				X	
14000/34	Objective semi plan 4x				X	
14000/35	Objective semi plan 10x				X	

Microscopes

109

Human Product Catalog 2022 | 2023



Accessories for LED Microscopes

Microscopes

110

REF	Description	Fluo	Premium	Advanced	Classic	Light
14000/36	Objective semi plan 20x				X	
14000/37	Objective semi plan 40x				X	
14000/38	Objective semi plan 100x				X	
14000/39	Objective plan 4x			X	X	
14000/40	Objective plan 10x			X	X	
14000/41	Objective plan 20x			X	X	
14000/42	Objective plan 40x			X	X	
14000/43	Objective plan 60x			X	X	
14000/44	Objective plan 100x			X	X	
14000/45	Objective plan infinity 4x	X	X			
14000/46	Objective plan infinity 10x	X	X			
14000/47	Objective plan infinity 20x	X	X			
14000/48	Objective plan infinity 40x		X			
14000/140	Objective plan infinity 60x	X	X			
14000/49	Objective plan infinity 100x	X	X			
14000/75	Oil	X	X	X	X	X
14000/82	Rubber eyeguard (set of 2)		X	X	X	X
14000/59	Simple achromatic phase contrast kit			X	X	X
14000/58	Simple polarisation kit		X	X	X	
14000/57	Simple polarisation kit					X
14000/68	Stage micrometer 1 mm / 100 division	X	X	X	X	X
14000/83	Transport box (plastic)			X	X	
14000/84	Transport box (plastic, light)					X
14000/13	Trinocular head			X	X	
14000/15	Trinocular infinity head		X			
14000/60	Turret achromatic phase contrast kit			X	X	
14000/62	Turret plan infinity phase contrast kit	X	X			
14000/61	Turret plan phase contrast kit			X	X	
14000/81	Yellow filter	X	X	X	X	X

Digital cameras for HumaScope Line (distributed by HUMAN)

REF	Description	Fluo	Premium	Advanced	Classic	Light
14000/866	Digital camera (USB 2) 5 megapixel	X	X	X	X	X
14000/867	Digital camera (USB 3) 5 megapixel	X	X	X	X	X
14000/877	Digital camera (USB 2) 10 megapixel	X	X	X	X	X
14000/878	Digital camera (USB 3) 10 megapixel	X	X	X	X	X

For operation with a trinocular, a C-mount adapter is necessary.