

AMH Control HumaCLIA SR

REF 82850

LOT 23001

Version 01 / 06 - 2023

Target Values / Sollwerte / Valores asignados / Valeurs cibles

🕒 2024-09-25

Link for latest version of target value sheet:

<https://www.human.de/clia/target-value-sheets/>



REF	Assay Parameter Ensayo	Valid from Gültig ab Válido a partir de	Units Einheiten Unidades	Target Sollwert Valor asignado Cible	SD SD SD	Range Bereich Intervalo	Target Sollwert Valor asignado Cible	SD SD SD	Range Bereich Intervalo
	Essai	Valable à partir de RGT Lot	Unité	ET	Intervalle	ET	Intervalle		
				Level 1			Level 2		
82150	AMH	22002	ng/ml	1,34	0,18	0,81 - 1,87	6,53	0,86	3,96 - 9,10

