

# HUMAN

Target Values

## BioRad Lyphocheck Immunoassay Plus Control

370 Trilevel 12x5ml  
371 Level 1 12x5ml  
372 Level 2 12x5ml  
373 Level 3 12x5ml  
370X MiniPak 3x5ml

Version 01/ 07-2023

REF

LOT 40420

Level 1 40421  
Level 2 40422  
Level 3 40423

IVD

2025-03-31

Link for latest version of target value sheet:

<https://www.human.de/cia/target-value-sheets/>



REF	Assay Parameter Ensayo Essai	Valid from Gültig ab Válido a partir de Valable à partir de RGT Lot	Units Einheiten Unidades Unité	Target	Range	Target	Range	Target	Range		
				Sollwert Valor asignado Cible	Bereich Intervalo Intervalle	Sollwert Valor asignado Cible	Bereich Intervalo Intervalle	Sollwert Valor asignado Cible	Bereich Intervalo Intervalle		
				Level 1 - 40421			Level 2 - 40422			Level 3 - 40423	
84600	TSH	20001	µIU/ml	0.47	0.31 - 0.63	6.32	4.11 - 8.53	41.4	26.9 - 55.8		
84620	fT3	20001	pg/ml	2.81	1.83 - 3.79	6.79	4.41 - 9.17	10.4	6.75 - 14.0		
84615	T3	20001	ng/ml	0.90	0.59 - 1.22	2.18	1.42 - 2.94	3.56	2.31 - 4.81		
84610	fT4	20001	pg/ml	7.99	5.19 - 10.80	19.9	12.9 - 26.8	38.6	25.1 - 52.1		
84605	T4	20001	ng/ml	54.8	35.6 - 73.9	115	75.1 - 156	157	102 - 212		
82120	FSH	20001	mIU/ml	7.71	5.01 - 10.40	23.8	15.5 - 32.2	53.6	34.9 - 72.4		
82115	LH	20001	mIU/ml	1.75	1.14 - 2.36	26.1	16.9 - 35.2	104	67.5 - 140		
83130	PSA	20001	ng/ml	0.95	0.62 - 1.28	5.42	3.52 - 7.32	36.0	23.4 - 48.6		
83125	fPSA	20001	ng/ml	0.48	0.31 - 0.65	2.28	1.48 - 3.08	11.3	7.36 - 15.3		
83100	AFP	20001	ng/ml	18.3	11.0 - 25.6	106	63.3 - 148	229	137 - 320		
82100	HCG	20001	mIU/ml	5.14	3.34 - 6.94	24.7	16.0 - 33.3	208	135 - 280		
82105	Progesterone	20001	ng/ml	0.19	0.12 - 0.26	11.1	7.19 - 14.90	29.1	18.9 - 39.3		
82110	Testosterone	20001	ng/ml	0.79	0.47 - 1.11	3.45	2.07 - 4.83	7.29	4.37 - 10.2		
85000	25-OH Vitamin D	22001	ng/ml	19.8	12.9 - 26.7	27.5	17.9 - 37.2	22.2	14.4 - 29.9		
83105	CEA	22001	ng/ml	2.95	1.92 - 3.98	17.3	11.2 - 23.3	33.6	21.8 - 45.3		
82130	Estradiol	22002	pg/ml	66.5	43.2 - 89.7	251	163 - 338	763	496 - 1029		
85049	Ferritin*	23001	ng/ml	74.0	48.1 - 99.9	192	125 - 259	458	298 - 619		
82125	Prolactin	22001	ng/ml	9.92	6.45 - 13.40	26.2	17.0 - 35.4	60.6	39.4 - 81.8		
85055	Vitamin B12	22002	pg/ml	442	287 - 597	717	466 - 968	1492	970 - 2015		

\* No CE - only for use outside Europe