


# PCT Control HumaCLIA SR

**REF** 85820

**LOT** 23001

Version 01 / 07 - 2023

Target Values / Sollwerte / Valores asignados / Valeurs cibles

 2026-10-12

Link for latest version of target value sheet:

<https://www.human.de/clia/target-value-sheets/>



<b>REF</b>	Assay Parameter Ensayo	Valid from Gültig ab Válido a partir de	Units Einheiten Unidades	Target Sollwert Valor asignado Cible	SD SD SD	Range Bereich Intervalo	Target Sollwert Valor asignado Cible	SD SD SD	Range Bereich Intervalo
	Essai	Valable à partir de <b>RGT</b> Lot	Unité		ET	Intervalle		ET	Intervalle
				Level 1			Level 2		
85020	PCT	23001	ng/ml	0.59	0.06	0.41 - 0.77	12.54	1.25	8.78 - 16.30