

# HbA1c liquidirect Calibrator Set

**LOT**

21002

Version 1 / 10 - 2021

**REF**

10776



29.02.2024

**Content** 4 x 0.5 ml

**Calibration Values / Kalibrationswerte / Valores de calibración / Valeurs de calibration**

			<b>CAL0</b>	<b>CAL1</b>	<b>CAL2</b>	<b>CAL3</b>	<b>CAL4</b>
Device Gerät Instrumento Equipement	Reference material/method Referenzmaterial/-methode Material/método de referencia Matériel/méthode de référence	Unit Einheit Unidad Unité	Assigned Value Zugewiesener Wert Valor asignado Valeur assignée	Assigned Value Zugewiesener Wert Valor asignado Valeur assignée	Assigned Value Zugewiesener Wert Valor asignado Valeur assignée	Assigned Value Zugewiesener Wert Valor asignado Valeur assignée	Assigned Value Zugewiesener Wert Valor asignado Valeur assignée
HumaStar 300SR	IFCC HbA1c system	mmol/mol Hb	0	29.0	73.0	115	140
		%	0	4.8	8.7	12.5	14.8
HumaStar 600	IFCC HbA1c system	mmol/mol Hb	17.0	41.0	81.0	126	159
		%	3.68	5.8	9.4	13.5	16.5
other	IFCC HbA1c system	mmol/mol Hb	0	33.0	77.0	119	147
		%	0	5.1	9.1	12.9	15.4