

# Standard RHEUMATOID FACTORS

Version 1/09 - 2022

REF

11361

LOT

22005



2024 - 02 - 29

Link for latest version of target value sheet:

<https://www.human.de/target-value-sheets/>

## Calibration Values / Kalibrationswerte / Valores de calibración / Valeurs de calibration

Device Gerät Instrumento Equipement	Reference material/method Referenzmaterial/-methode Material/método de referencia Matériel/méthode de référence	Unit Einheit Unidad Unité	Assigned Value Zugewiesener Wert Valor asignado Valeur assignée
<b>HUMAN RHEUMATOID FACTORS - Immunoturbidimetry Latex enhanced</b>			
Multipurpose reagent on HumaStar 100 & 200	Adjusted to 1 <sup>st</sup> British Standard NIBSC 0064/2	IU/ml	114
		klU/l	114
Multipurpose reagent on HumaLyzer 4000	Adjusted to 1 <sup>st</sup> British Standard NIBSC 0064/2	IU/ml	114
		klU/l	114
Multipurpose reagent for other	Adjusted to 1 <sup>st</sup> British Standard NIBSC 0064/2	IU/ml	114
		klU/l	114
HumaStar System reagent on HumaStar 300SR	Adjusted to 1 <sup>st</sup> British Standard NIBSC 0064/2	IU/ml	114
		klU/l	114
HumaStar System reagent on HumaStar 600	Adjusted to 1 <sup>st</sup> British Standard NIBSC 0064/2	IU/ml	103
		klU/l	103

CA\_RF\_11361

Rev. 001 | valid of 17.03.2021

INF 1136102 - Version 1/09 - 2022



**Human**

Human Gesellschaft für Biochemica und Diagnostica mbH

Max-Planck-Ring 21 · 65205 Wiesbaden · Germany

Telefon +49 6122-9988-0 · Telefax +49 612-9988-100 · e-Mail human@human.de