

Standard RHEUMATOID FACTORS

Version 1/02 - 2023

REF

11361

LOT

23004



2024 - 11 - 30

Calibration Values / Kalibrationswerte / Valores de calibración / Valeurs de calibration

Device Gerät Instrumento Equipement	Reference material/method Referenzmaterial/-methode Material/método de referencia Matériel/méthode de référence	Unit Einheit Unidad Unité	Assigned Value Zugewiesener Wert Valor asignado Valeur assignée
RHEUMATOID FACTORS			
Multipurpose reagent, REF 11261PA			
HumaStar 100 & 200	Adjusted to 1 st British Standard NIBSC 0064/2	IU/ml	114
		kIU/l	114
HumaLyzer 4000	Adjusted to 1 st British Standard NIBSC 0064/2	IU/ml	114
		kIU/l	114
Other	Adjusted to 1 st British Standard NIBSC 0064/2	IU/ml	114
		kIU/l	114
System reagents, REF 11261300, 11261600			
HumaStar 300SR	Adjusted to 1 st British Standard NIBSC 0064/2	IU/ml	114
		kIU/l	114
HumaStar 600	Adjusted to 1 st British Standard NIBSC 0064/2	IU/ml	103
		kIU/l	103

CA_RF_11361

Rev. 001 | valid of 17.03.2021

INF 1136102 - Version 1/02 - 2023

Human

Human Gesellschaft für Biochemica und Diagnostica mbH
Max-Planck-Ring 21 · 65205 Wiesbaden · Germany
Telefon +49 6122-9988-0 · Telefax +49 612-9988-100 · e-Mail human@human.de