

Standard RHEUMATOID FACTORS

Version 1/11 - 2023

REF

11361

LOT

24001



2025 - 04 - 30

Link for latest version of target value sheet:

<https://www.human.de/target-value-sheets/>

Calibration Values / Kalibrationswerte / Valores de calibración / Valeurs de calibration

Device Gerät Instrumento Equiptement	Reference material/method Referenzmaterial/-methode Material/método de referencia Matériel/méthode de référence	Unit Einheit Unidad Unité	Assigned Value Zugewiesener Wert Valor asignado Valeur assignée
RHEUMATOID FACTORS			
Multipurpose reagent, REF 11261PA			
HumaStar 100 & 200 LOT ≥ 23501	In house Master calibrator traceable to 1st British Standard NIBSC 0064/2	IU/ml	180
		kIU/l	180
HumaLyzer 4000 LOT ≥ 23501	In house Master calibrator traceable to 1st British Standard NIBSC 0064/2	IU/ml	180
		kIU/l	180
Other LOT ≥ 23501	In house Master calibrator traceable to 1st British Standard NIBSC 0064/2	IU/ml	180
		kIU/l	180
System reagents, REF 11261300, 11261600			
HumaStar 300SR LOT ≥ 23501	In house Master calibrator traceable to 1st British Standard NIBSC 0064/2	IU/ml	180
		kIU/l	180
HumaStar 600 LOT ≥ 23501	In house Master calibrator traceable to 1st British Standard NIBSC 0064/2	IU/ml	180
		kIU/l	180

