


HbA1c liquidirect Controls

REF 10775

Content 2 x 0.5 ml

LOT 21001
 29.02.2024

Version 1 / 26-07-2021

Target Values / Sollwerte / Valores Asignados/ Valeurs Cibles

Unit
Einheit
Unidad
Unité

Control Blood normal CBN	
Mean	Range
Mittelwert	Bereich
Media	Rango
Moyenne	Plage de valeurs

Control Blood abnormal CBA	
Mean	Range
Mittelwert	Bereich
Media	Rango
Moyenne	Plage de valeurs

Device Gerät Instrumento Equipement	HumaStar 300SR	mmol/mol Hb
		%
	HumaStar 600	mmol/mol Hb
		%
other	mmol/mol Hb	
	%	

38,8	34,0	—	43,5
5,64	5,21	—	6,07
46,5	45,7	—	47,3
6,34	6,26	—	6,41
39,3	32,8	—	45,8
5,69	5,10	—	6,27

118	112	—	125
12,8	12,2	—	13,4
124	121	—	128
13,3	13,0	—	13,7
118	112	—	125
12,8	12,2	—	13,4