

# HbA1c liquidirect Controls

REF 10775

Content 2 x 0.5 ml

LOT 22001  
 29.02.2024

Version 1 / 26-07-2021

## Target Values / Sollwerte / Valores Asignados/ Valeurs Cibles

|         |
|---------|
| Unit    |
| Einheit |
| Unidad  |
| Unité   |

| Control Blood normal<br>[CBN] |                  |
|-------------------------------|------------------|
| Mean                          | Range            |
| Mittelwert                    | Bereich          |
| Media                         | Rango            |
| Moyenne                       | Plage de valeurs |

| Control Blood abnormal<br>[CBA] |                  |
|---------------------------------|------------------|
| Mean                            | Range            |
| Mittelwert                      | Bereich          |
| Media                           | Rango            |
| Moyenne                         | Plage de valeurs |

| Device<br>Gerät<br>Instrumento<br>Equipement | HumaStar 300SR | mmol/mol Hb |
|--|----------------|-------------|
|  |                | %           |
|  | HumaStar 600   | mmol/mol Hb |
|  |                | %           |
|  | other          | mmol/mol Hb |
|  |                | %           |

|      |      |   |      |
|------|------|---|------|
| 38.8 | 34.0 | — | 43.5 |
| 5.64 | 5.21 | — | 6.07 |
| 46.5 | 45.7 | — | 47.3 |
| 6.34 | 6.26 | — | 6.41 |
| 39.3 | 32.8 | — | 45.8 |
| 5.69 | 5.10 | — | 6.27 |

|      |      |   |      |
|------|------|---|------|
| 118  | 112  | — | 125  |
| 12.8 | 12.2 | — | 13.4 |
| 124  | 121  | — | 128  |
| 13.3 | 13.0 | — | 13.7 |
| 118  | 112  | — | 125  |
| 12.8 | 12.2 | — | 13.4 |