


HbA1c liquidirect Controls

LOT 23002

Version 1 /12-06 - 2023

REF 10775

 2025-04-30

Link for latest version of target value sheet:

<https://www.human.de/target-value-sheets/>

Target Values / Sollwerte / Valores Asignados/ Valeurs Cibles

Device Gerät Instrumento Equipement		Unit Einheit Unidad Unité	Control Blood normal CBN				Control Blood abnormal CBA			
			Mean Mittelwert Media Moyenne	Range Bereich Rango Plage de valeurs			Mean Mittelwert Media Moyenne	Range Bereich Rango Plage de valeurs		
HumaStar 300SR	mmol/mol Hb		36.1	29.6	—	42.6	108	101	—	114
	%		5.5	4.9	—	6.1	12.0	11.4	—	12.6
HumaStar 600	mmol/mol Hb		44.4	40.3	—	48.5	120	113	—	128
	%		6.2	5.8	—	6.6	13.1	12.5	—	13.9
other	mmol/mol Hb		36.5	29.4	—	43.6	112	97.3	—	126
	%		5.5	4.8	—	6.1	12.4	11.1	—	13.7

CS-HbA1c

Rev. 006 | valid of 31.03.2022

INF 1077502 - Version 1 /12-06 - 2023



Human Gesellschaft für Biochemica und Diagnostica mbH

Max-Planck-Ring 21 · 65205 Wiesbaden · Germany

Telefon +49 6122-9988-0 · Telefax +49 612-9988-100 · e-Mail human@human.de