

HUMAN

Target Values

BioRad Lyphocheck Tumor Marker Plus Control

REF 368X, **LOT** 54690 Level 1 - 3  2023-06-30

Version 2

REF		Units	Target	Range	Target	Range	Target	Range
		Einheiten	Sollwert	Bereich	Sollwert	Bereich	Sollwert	Bereich
		Unités	Valor meta	Rango	Valor meta	Rango	Valor meta	Rango
		Unités	Valeur souhaitée	Marge	Valeur souhaitée	Marge	Valeur souhaitée	Marge
HUMAN Tumor Markers								
			Level 1 - 54691		Level 2 - 54692		Level 3 - 54693	
52010	AFP	ng/ml	8.79	5.89 - 11.7	84.7	56.7 - 113	264	177 - 351
52020	CEA	ng/ml	2.94	1.97 - 3.91	26.8	18.0 - 35.7	91.9	61.6 - 122
52050	CA125	U/ml	24.9	16.7 - 33.1	71.4	47.9 - 95.0	221	148 - 294
52080	CA15-3	U/ml	27.2	18.2 - 36.2	53.6	35.9 - 71.3	92.6	62.0 - 123
52060	CA19-9	U/ml	16.7	11.2 - 22.3	33.3	22.3 - 44.2	81.7	54.7 - 109
52035	fPSA	ng/ml	<0.5	-	2.72	1.83 - 3.62	>10	-
52030	PSA	ng/ml	0.10	0.06 - 0.15	6.07	3.95 - 8.2	23.9	15.5 - 32.2



#54690
09-2018-02

Human Gesellschaft für Biochemica und Diagnostica mbH
Max-Planck-Ring 21 - D-65205 Wiesbaden - Germany
Tel.: +496122 99880 Fax: +496122 9988100 - e-mail: human@human.de