

# HUMAN

Target Values

## BioRad Lyphocheck Tumor Marker Plus Control

**REF** 368X, **LOT** 74610 Level 1 - 3  2024-09-30

REF		Units	Target	Range	Target	Range	Target	Range
		Einheiten Unitades Unités	Sollwert Valor meta Valeur souhaitée	Bereich Rango Marge	Sollwert Valor meta Valeur souhaitée	Bereich Rango Marge	Sollwert Valor meta Valeur souhaitée	Bereich Rango Marge
<b>HUMAN Tumor Markers</b>								
			Level 1 - 74611		Level 2 - 74612		Level 3 - 74613	
52010	AFP	ng/ml	8.37	5.61 - 11.1	91.6	61.4 - 122	257	172 - 342
52020	CEA	ng/ml	2.68	1.80 - 3.57	23.5	15.8 - 31.3	85.1	57.0 - 113
52050	CA125	U/ml	18.9	11.6 - 26.3	59.9	40.1 - 79.6	212	142 - 282
52080	CA15-3	U/ml	35.2	23.6 - 46.8	79.4	53.2 - 105.6	163	109 - 217
52060	CA19-9	U/ml	13.9	9.29 - 18.4	24.3	16.3 - 32.3	80.1	53.7 - 107
52035	fPSA	ng/ml	<0,3		2.71	1.82 - 3.60	>10	
53040	PSA	ng/ml	0.21	0.11 - 0.29	7.97	5.18 - 10.7	27.9	18.1 - 37.6



Lypho 74610  
09-2018-02

Human Gesellschaft für Biochemica und Diagnostica mbH  
Max-Planck-Ring 21 - D-65205 Wiesbaden - Germany  
Tel.: +496122 99880 Fax: +496122 9988100 - e-mail: human@human.de