

HUMAN

Target Values

BioRad Lyphocheck Tumor Marker Plus Control

Version 01/11 - 2021

367 Level 1 6x2ml Level 1 54691
 [REF] 368 Level 2 6x2ml [LOT] 54690 Level 2 54692 [IVD] 2023-06-30
 369 Level 3 6x2ml Level 3 54693
 368X Trilevel MiniPak 3x2ml

Link for latest version of target value sheet:
<https://www.human.de/clia/target-value-sheets/>



[REF]	Assay Parameter Ensayo	Valid from Gültig ab Válido a partir de RGT Lot	Units Einheiten Unidades Unité	Target Sollwert Valor asignado Cible	Range Bereich Intervalo Intervalle	Target Sollwert Valor asignado Cible	Range Bereich Intervalo Intervalle	Target Sollwert Valor asignado Cible	Range Bereich Intervalo Intervalle
					Level 1 - 54691	Level 2 - 54692		Level 3 - 54693	
83130	PSA	20001	ng/ml	0.07	0.01 - 0.13	3.22	2.12 - 4.37	15.78	10.37 - 21.42
83125	fPSA	20001	ng/ml	0.03	0.00 - 0.06	2.75	1.81 - 3.73	12.66	8.32 - 17.18
83100	AFP	20001	ng/ml	9.66	6.47 - 12.85	74.04	49.61 - 98.47	207.49	139.02 - 275.96

#54690
 Rev. 01 | valid of 10.11.2021

