

HUMAN

Target Values

BioRad Lyphocheck Immunoassay Plus Control

370 Trilevel 12x5ml
 371 Level 1 12x5ml
 372 Level 2 12x5ml
 373 Level 3 12x5ml
 370X MiniPak 3x5ml

Level 1 40391
 Level 2 40392
 Level 3 40393

REF

LOT 40390

IVD

2022-12-31

Link for latest version of target value sheet:

<https://www.human.de/cia/target-value-sheets/>



REF	Assay Parameter Ensayo	Valid from Gültig ab Válido a partir de RGT Lot	Units Einheiten Unidades Unité	Target Sollwert Valor asignado Cible	Range Bereich Intervalo Intervalle	Target Sollwert Valor asignado Cible	Range Bereich Intervalo Intervalle	Target Sollwert Valor asignado Cible	Range Bereich Intervalo Intervalle
				Level 1 - 40391	Level 2 - 40392		Level 3 - 40393		
84600	TSH	20001	µIU/ml	0,660	0,430 - 0,900	5,89	3,87 - 7,99	37,3	24,5 - 50,6
84620	FT3	20001	pg/ml	2,50	1,64 - 3,39	5,93	3,90 - 8,05	8,46	5,56 - 11,5
84610	FT4	20001	pg/ml	5,83	3,83 - 7,91	13,0	8,55 - 17,6	19,8	13,0 - 26,9
82120	FSH	20001	mIU/ml	5,50	3,62 - 7,46	22,8	15,0 - 30,9	51,8	34,1 - 70,3
84605	T4	20001	ng/ml	50,0	32,9 - 67,8	94,3	62,0 - 128	112	73,3 - 151
83130	PSA*	20001	ng/ml	1,06	0,700 - 1,44	4,54	2,98 - 6,16	27,0	17,7 - 36,6

*Conformity assessment pending