

HUMAN

Target Values

BioRad Lyphocheck Immunoassay Plus Control

370 Trilevel 12x5ml
371 Level 1 12x5ml
REF 372 Level 2 12x5ml
373 Level 3 12x5ml
370X MiniPak 3x5ml

Level 1 40401
Level 2 40402
Level 3 40403

Version 03/ 11 - 2021

LOT 40400

IVD

2023-10-31

Link for latest version of target value sheet:

<https://www.human.de/cia/target-value-sheets/>



REF	Assay Parameter Ensayo	Valid from Gülig ab Valable à partir de RST Lot	Units Einheiten Unidades	Target Sollwert Valor asignado Cible	Range Bereich Intervalo Intervalle	Target Sollwert Valor asignado Cible	Range Bereich Intervalo Intervalle	Target Sollwert Valor asignado Cible	Range Bereich Intervalo Intervalle
84600	TSH	20001	µIU/ml	0.68	0.45 - 0.92	6.08	4.00 - 8.25	35.14	23.10 - 47.69
84620	FT3	20001	pg/ml	2.49	1.64 - 3.38	6.46	4.25 - 8.77	9.59	6.30 - 13.02
84615	T3	20001	ng/ml	0.80	0.53 - 1.09	2.14	1.41 - 2.90	3.35	2.20 - 4.55
84610	FT4	20001	pg/ml	6.00	3.94 - 8.14	14.10	9.27 - 19.14	26.59	17.48 - 36.09
84605	T4	20001	ng/ml	53.66	35.27 - 72.82	112.6	73.99 - 152.8	154.29	101.42 - 209.39
82120	FSH	20001	mIU/ml	5.42	3.56 - 7.36	23.43	15.40 - 31.80	43.29	28.46 - 58.75
82115	LH	20001	mIU/ml	2.13	1.40 - 2.89	32.18	21.15 - 43.7	118.25	77.73 - 160.48
83130	PSA	20001	ng/ml	1.08	0.710 - 1.47	4.20	2.76 - 5.70	27.17	17.86 - 36.87
83125	fPSA	20001	ng/ml	0.36	0.24 - 0.49	1.96	1.29 - 2.66	11.85	7.79 - 16.08
83100	AFP	20001	ng/ml	16.61	9.97 - 23.25	87.8	52.68 - 122.92	191.93	115.16 - 268.70
82100	HCG	20001	mIU/ml	5.13	3.37 - 6.96	23.16	15.22 - 31.43	161.89	106.41 - 219.71
82105	Progesterone*	20001	ng/ml	<0,1	<0,1 - <0,1	12.69	8.34 - 17.22	33.04	21.72 - 44.84
82110	Testosterone*	20001	ng/ml	0.90	0.54 - 1.26	3.41	2.05 - 4.77	6.09	4.14 - 9.66

*Conformity assessment pending

#40400
Rev. 03 | valid of 10.11.2021