


Battery system and solar panel for all HUMAN devices – at least 2 hours battery operation



Clinical Chemistry	Hours
HumaLyzer 4000 / Primus	> 8
HumaStar 100 / 200 / 300SR*	> 4
HumaLyte Plus3/5	> 12



Hematology	Hours
HumaCount 5D ^{CRP} / HumaCount 5D	> 5
HumaCount 5L	> 3
HumaCount 80 ^{TS} / 30 ^{TS}	> 12
HumaSRate 24 ^{PT}	> 12



Clinical Chemistry	Hours
HumaReader HS / Single plus	> 12
Elisys Uno*	> 4
Elisys Duo	> 3
Combiwash	> 5



Hemostasis	Hours
HumaClot Junior	> 48
HumaClot Duo plus / Quattro	> 8
HumaClot Pro	> 4

Urinalysis	Hours
Combilyzer 13	> 8



Laboratory Line	Hours
HumaPure	> 1000
HuMax Micro / HumaRock / HumaRoll	> 48
HumaTherm	> 8
HumaCube / HumaTwist	> 6
HuMax 3K	> 5
HuMax 4K / 5K / HCT	> 2

Microscopes	Hours
HumaScope Light / Classic / Advanced	> 1000

Diabetes DX	Hours
HumaMeter A1c	> 24
HumaNex A1c Variant	> 4

Autoimmune DX	Hours
HumaBlot 44 ^{FA}	> 5



Molecular DX	Hours
HumaLoop T / M	> 4
HumaTurb C+A	> 24
HumaHeat	< 4

* Analyzer without PC