

HumaType REAGENTS

| Reagents for Blood Grouping
& Compatibility Testing

Immunoematology

- > Strong and robust agglutinations
- > Meets high European Quality Standards CE

Human

Diagnostics Worldwide

HumaType REAGENTS

for Blood Grouping & Compatibility Testing



AVAILABLE AT HUMAN

Line of potent blood grouping reagents with high antigen affinity. Reliable reagents for compatibility testing and antibody screening. Suited for traditional techniques. High quality, CE-marked products.

ABO REAGENTS

- > Highly avid monoclonal reagents
- > Excellent sensitivity – even for the weak phenotypes A₃, A_x and B_x
- > Superior specificity – no false reactions with polyagglutinable cells
- > Suited for the tube, slide and microplate technique

		Cat. No.
HumaType ANTI-A	3 x 10 ml	71012
HumaType ANTI-B	3 x 10 ml	71022
HumaType ANTI-A,B	3 x 10 ml	71032



RHESUS REAGENTS

- > Highly avid monoclonal reagents
- > Suited for the tube, slide and microplate technique

HumaType ANTI-D (IgM)

- > Detects all D antigens – except category D^v cells
- > The ideal reagent for typing patient specimens

3 x 10 ml Cat. No. 71042

HumaType ANTI-D (IgM/IgG)

- > Blend of IgM & IgG anti-D – suitable for patient and donor typing
- > Detects all D antigens – in direct or indirect agglutination tests

3 x 10 ml Cat. No. 71052

COMPATIBILITY TESTING/ANTIBODY DETECTION

HumaType ANTI-HUMAN GLOBULIN

- > Optimized blend of rabbit anti-IgG, anti-C3b / -C3d
- > For direct and indirect antiglobulin tests
- > Dyed green for easy visualization

3 x 10 ml Cat. No. 71062

HumaType ALBUMIN 22%

- > Potentiator to enhance antigen/antibody reactions
- > Polymer enhanced bovine serum albumin

3 x 10 ml Cat. No. 71072

Human

Diagnostics Worldwide