

NEW

HumaTube Line

HUMAN's Quality
Sterile Blood Collection Devices

- > Reliable
- > Safe
- > Convenient

BLOOD COLLECTION



Human

Diagnostics Worldwide

HumaTube Line

Benefits for blood collection and the clinical laboratory

PET Vacuum Tubes

Unbreakable, glass-clear PET

- > Safe
- > Light weight: Lower transport and disposal costs

Colour-coded caps



Evacuated and sterile

- > Accurate filling for reliable results
- > Convenient use

- > Clear identification
- > Safe sealing (cap with inner stopper part)
- > Internationally accepted colour code

STERILE

Aluminium bag for citrate tubes

- > Extended shelf life



Serum Tubes

REF

No Additive (Z)	5 ml	12 x 100 tubes	73005/12
	7 ml	12 x 100 tubes	73010/12
	10 ml	8 x 100 tubes	73013/08 NEW
Serum C/A (Z)	5 ml	12 x 100 tubes	73023/12 NEW
	7 ml	12 x 100 tubes	73020/12
Serum Gel-C/A (Z)	5 ml	12 x 100 tubes	73030/12

Plasma Tubes

Li-Heparin (LH)	3 ml	12 x 100 tubes	73050/12
Na-Citrate 3.2 % (9NC)	2.7 ml	48 x 25 tubes	73060/12
Glucose (FH)	2 ml	12 x 100 tubes	73070/12

Whole-Blood Tubes

K ₃ -EDTA (K ₃ E)	3 ml	12 x 100 tubes	73040/12
	7 ml	12 x 100 tubes	73045/12 NEW

ESR see HumaSed Line

Needles and Accessories

NEW



Latex free and sterile

- > patient safety

Sharp & smooth edges

- > painless penetration

STERILE

Colour coded

- > clear identification of lumen size
- > green (21G)
- > black (22G)
- > blue (23G)

Double wings

- > safe puncture

Holder

- > easy connection to rubber stoppers

Winged blood collection set

REF

HumaTube Wing Set 21G, 3/4"	0.80 x 20 mm	40 x 50 pcs	73111
HumaTube Wing Set 22G, 3/4"	0.70 x 20 mm	40 x 50 pcs	73112
HumaTube Wing Set 23G, 3/4"	0.60 x 20 mm	40 x 50 pcs	73113

Pentype blood collection needle

HumaTube Pen 21G, 1 1/2"	0.80 x 38 mm	40 x 100 pcs	73131
HumaTube Pen 22G, 1 1/2"	0.70 x 38 mm	40 x 100 pcs	73132
HumaTube Pen 23G, 1 1/2"	0.60 x 38 mm	40 x 100 pcs	73133

Blood collection tube holder

HumaTube Holder	2000 pcs	73100
-----------------	----------	-------