


# Diagnosis of APS

## In Case of Suspected Antiphospholipid Syndrome

### AUTOIMMUNE DIAGNOSTICS



# ANTI-PHOSPHOLIPID ANTIBODIES

## ELISA



New

**PHOSPHOLIPID-ANTIBODIES SCREEN** 96 tests ITC59070 REF

ELISA for the Quantitative Determination of Anti-Phospholipid Antibodies (IgG and IgM) Cardiolipin,  $\beta_2$ GP1, Phosphatidylserine, Phosphatidylethanolamine, Phosphatidic acid, phosphatidylinositol

**CARDIOLIPIN-ANTIBODIES SCREEN** 96 tests ITC59076

ELISA for the Quantitative Determination of Anti-Cardiolipin Antibodies (Ig(GAM)) Cardiolipin + cofactor  $\beta_2$ GP1

**CARDIOLIPIN-ANTIBODIES COMBI** 96 tests ITC59082

ELISA for the Quantitative Determination of Anti-Cardiolipin Antibodies (IgG and IgM) Cardiolipin + cofactor  $\beta_2$ GP1

**CARDIOLIPIN-ANTIBODIES IgG** 96 tests ITC59071

ELISA for the Quantitative Determination of Anti-Cardiolipin Antibodies (IgG) Cardiolipin + cofactor  $\beta_2$ GP1

**CARDIOLIPIN-ANTIBODIES IgM** 96 tests ITC59081

ELISA for the Quantitative Determination of Anti-Cardiolipin Antibodies (IgM) Cardiolipin + cofactor  $\beta_2$ GP1

**PHOSPHATIDYLSERINE-ANTIBODIES SCREEN** 96 tests ITC59027

ELISA for the Quantitative Determination of Anti-Phosphatidylserine Antibodies (Ig(GAM)) Phosphatidylserine + cofactor  $\beta_2$ GP1

**PHOSPHATIDYLSERINE-ANTIBODIES IgG** 96 tests ITC59011

ELISA for the Quantitative Determination of Anti-Phosphatidylserine Antibodies (IgG) Phosphatidylserine + cofactor  $\beta_2$ GP1

**PHOSPHATIDYLSERINE-ANTIBODIES IgM** 96 tests ITC59021 REF

ELISA for the Quantitative Determination of Anti-Phosphatidylserine Antibodies (IgM) Phosphatidylserine + cofactor  $\beta_2$ GP1

**PHOSPHATIDYLETHANOLAMINE-ANTIBODIES SCREEN**

ELISA for the Quantitative Determination of Anti-Phosphatidylethanolamine Antibodies (Ig(GAM)) Phosphatidylethanolamine

**$\beta_2$ -GLYCOPROTEIN I-ANTIBODIES SCREEN** 96 tests ITC59050

ELISA for the Quantitative Determination of Anti- $\beta_2$ -Glycoprotein-I Antibodies (Ig(GAM))  $\beta_2$ GP1

**$\beta_2$ -GLYCOPROTEIN I-ANTIBODIES IgG** 96 tests ITC59150

ELISA for the Quantitative Determination of Anti- $\beta_2$ -Glycoprotein I Antibodies (IgG)  $\beta_2$ GP1

**$\beta_2$ -GLYCOPROTEIN I-ANTIBODIES IgM** 96 tests ITC59250

ELISA for the Quantitative Determination of Anti- $\beta_2$ -Glycoprotein I Antibodies (IgM)  $\beta_2$ GP1

**PROTHROMBIN-ANTIBODIES SCREEN** 96 tests ITC59450

ELISA for the Quantitative Determination of Anti-Prothrombin Antibodies (Ig(GAM)) Prothrombin (factor II), phosphatidylserine

**ANNEXIN V-ANTIBODIES SCREEN** 96 tests ITC59550

ELISA for the Quantitative Determination of Anti-Annexin V Antibodies (Ig(GAM)) Annexin V

