

**HUMAN**  
Target Values  
**BioRad Lyphocheck Tumor Marker Plus Control**  
**REF 368X, LOT 54620 Level 1 - 3**

REF		Units	Target	Range	Target	Range	Target	Range
		Einheiten Unitades Unités	Sollwert Valor meta Valeur souhaitée	Bereich Rango Marge	Sollwert Valor meta Valeur souhaitée	Bereich Rango Marge	Sollwert Valor meta Valeur souhaitée	Bereich Rango Marge
<b>HUMAN Tumor Markers</b>								
			Level 1 - 54621		Level 2 - 54622		Level 3 - 54623	
52010	AFP	ng/ml	6.9	4.6 - 9.2	65.29	43.8 - 86.8	203	136 - 270
52020	CEA	ng/ml	4.02	2.69 - 5.34	22.32	15.0 - 29.7	81.5	54.6 - 108
52050	CA125	U/ml	21.69	14.5 - 28.9	55.4	37.1 - 73.7	165.2	111 - 220
52080	CA15-3	U/ml	23.96	16 - 31.9	58.2	39.0 - 77	114.1	76.5 - 152
52060	CA19-9	U/ml	17.1	11.4 - 22.7	35.7	23.9 - 47.5	108	72.4 - 144
52035	fPSA	ng/ml	<0.5	-	2.98	2.00 - 3.96	>10	-
53040	PSA	ng/ml	0.23	0.1 - 0.33	9.7	8.2 - 11.0	38	34.0 - 42



Lypho 54620  
04-2016-1

Human Gesellschaft für Biochemica und Diagnostica mbH  
Max-Planck-Ring 21 - D-65205 Wiesbaden - Germany  
Tel.: +496122 99880 Fax: +496122 9988100 - e-mail: human@human.de