

HUMAN

Target Values

BioRad Lyphocheck Tumor Marker Plus Control

[REF] 368X, [LOT] 54630 Level 1 - 3 ☒ 2019-09-30

[REF]		Units Einheiten Unitades Unitées	Target Sollwert Valor meta Valeur souhaitée	Range Bereich Rango Marge	Target Sollwert Valor meta Valeur souhaitée	Range Bereich Rango Marge	Target Sollwert Valor meta Valeur souhaitée	Range Bereich Rango Marge
HUMAN Tumor Markers								
			Level 1 - 54621		Level 2 - 54622		Level 3 - 54623	
52010	AFP	ng/ml	7.14	4.78 - 9.49	66.7	44.7 - 88.7	202	135 - 268
52020	CEA	ng/ml	3.66	2.45 - 4.86	26.0	17.4 - 34.6	98.8	66.2 - 131
52050	CA125	U/ml	18.7	12.5 - 24.9	59.6	39.9 - 79.3	131	88.0 - 175
52080	CA15-3	U/ml	24.1	16.2 - 32.1	54.3	36.4 - 72.3	96.0	64.3 - 128
52060	CA19-9	U/ml	16.8	11.2 - 22.3	38.2	25.6 - 50.8	91.5	61.3 - 122
52035	fPSA	ng/ml	<0.5	-	1.92	1.29 - 2.55	8.99	6.03 - <10
53040	PSA	ng/ml	0.08	0.00 - 0.12	5.88	5.10 - 6.66	27.7	23.5 - 31.9



Lypho 54630
08-2016-1

Human Gesellschaft für Biochemica und Diagnostica mbH
Max-Planck-Ring 21 - D-65205 Wiesbaden - Germany
Tel.: +496122 99880 Fax: +496122 9988100 - e-mail: human@human.de