

# HUMAN

Target Values

## BioRad Lyphocheck Tumor Marker Plus Control

**REF** 368X, **LOT** 54650 Level 1 - 3  2020-07-31

REF		Units	Target	Range	Target	Range	Target	Range
		Einheiten	Sollwert	Bereich	Sollwert	Bereich	Sollwert	Bereich
		Unitades	Valor meta	Rango	Valor meta	Rango	Valor meta	Rango
		Unités	Valeur souhaitée	Marge	Valeur souhaitée	Marge	Valeur souhaitée	Marge
<b>HUMAN Tumor Markers</b>								
			Level 1 - 54641		Level 2 - 54642		Level 3 - 54643	
52010	AFP	ng/ml	9.55	6.40 - 12.7	77.1	51,7 - 103	234	157 - 312
52020	CEA	ng/ml	3.54	2.37 - 4.71	23.3	15.6 - 31.0	76.2	51,0 - 101
52050	CA125	U/ml	19.6	13.1 - 26.0	55.6	37.3 - 74.0	168	112 - 223
52080	CA15-3	U/ml	18.0	12.0 - 23.9	34.5	23.1 - 45.8	67.2	45.0 - 89.4
52060	CA19-9	U/ml	15.0	10.0 - 19.9	31.0	20.8 - 41.3	83.4	55.9 - 111
52035	fPSA	ng/ml	<0.5	-	2.82	1.89 - 3.76	>10	-
53040	PSA	ng/ml	0.11	0.06 - 0.16	7.53	4.90 - 10.2	33.9	22.0 - 45.7



Lypho 54650  
09-2017-01

Human Gesellschaft für Biochemica und Diagnostica mbH  
Max-Planck-Ring 21 - D-65205 Wiesbaden - Germany  
Tel.: +496122 99880 Fax: +496122 9988100 - e-mail: human@human.de