

# HUMAN

Target Values

## BioRad Lyphocheck Tumor Marker Plus Control

**REF** 368X, **LOT** 54670 Level 1 - 3  2022-03-31

Version 2

REF		Units	Target	Range	Target	Range	Target	Range
		Einheiten	Sollwert	Bereich	Sollwert	Bereich	Sollwert	Bereich
		Unités	Valor meta	Rango	Valor meta	Rango	Valor meta	Rango
		Unités	Valeur souhaitée	Marge	Valeur souhaitée	Marge	Valeur souhaitée	Marge
<b>HUMAN Tumor Markers</b>								
			Level 1 - 54671		Level 2 - 54672		Level 3 - 54673	
52010	AFP	ng/ml	7.89	5.29 - 10.5	87.7	58.7 - 117	240	161 - 320
52020	CEA	ng/ml	2.93	1.96 - 3.90	28.1	18.8 - 37.4	100.9	67.6 - 134
52050	CA125	U/ml	22.8	15.3 - 30.3	63.9	42.8 - 85.0	199	133 - 264
52080	CA15-3	U/ml	24.8	16.6 - 33.0	48.9	31.2 - 66.6	86.9	53.6 - 120
52060	CA19-9	U/ml	18.1	12.1 - 24.1	38.1	25.5 - 50.6	86.1	57.7 - 115
52035	fPSA	ng/ml	<0.5	-	3.18	2.13 - 4.23	>10	-
52030	PSA	ng/ml	0.21	0.12 - 0.31	7.42	4.82 - 10.0	28	18.2 - 37.9



Lypho 54660  
08-2018-01

Human Gesellschaft für Biochemica und Diagnostica mbH  
Max-Planck-Ring 21 - D-65205 Wiesbaden - Germany  
Tel.: +496122 99880 Fax: +496122 9988100 - e-mail: human@human.de