

# HUMAN

Target Values

## BioRad Lyphocheck Tumor Marker Plus Control

**REF** 368X, **LOT** 54680 Level 1 - 3  2022-08-31

Version 2

REF		Units	Target	Range	Target	Range	Target	Range
		Einheiten	Sollwert	Bereich	Sollwert	Bereich	Sollwert	Bereich
		Unités	Valor meta	Rango	Valor meta	Rango	Valor meta	Rango
		Unités	Valeur souhaitée	Marge	Valeur souhaitée	Marge	Valeur souhaitée	Marge
<b>HUMAN Tumor Markers</b>								
			Level 1 - 54681		Level 2 - 54682		Level 3 - 54683	
52010	AFP	ng/ml	7.83	5.25 - 10.4	81.0	54.3 - 108	233	153 - 310
52020	CEA	ng/ml	3.67	2.46 - 4.88	28.9	19.3 - 38.4	94.1	63.1 - 125
52050	CA125	U/ml	26.5	17.8 - 35.2	74.6	50.0 - 99.2	235	157 - 313
52080	CA15-3	U/ml	30.4	20.4 - 40.5	55.7	37.3 - 74.1	104.7	70.2 - 139
52060	CA19-9	U/ml	17.0	11.4 - 22.6	34.5	23.1 - 45.9	89.1	59.7 - 118
52035	fPSA	ng/ml	<0.5	-	3.67	2.46 - 4.88	>10	-
52030	PSA	ng/ml	0.06	0.03 - 0.08	7.47	4.86 - 10.1	26.5	17.2 - 35.8



#54680  
09-2018-02

Human Gesellschaft für Biochemica und Diagnostica mbH  
Max-Planck-Ring 21 - D-65205 Wiesbaden - Germany  
Tel.: +496122 99880 Fax: +496122 9988100 - e-mail: human@human.de