

# Immunoassay Multi Control HumaCLIA SR

REF

84850

LOT

22001

Version 02 / 03 - 2023

Target Values / Sollwerte / Valores asignados / Valeurs cibles



2025-06-24

Link for latest version of target value sheet:

<https://www.human.de/clia/target-value-sheets/>



REF	Assay Parameter Ensayo  Essai	Valid from Gültig ab Válido a partir de  Valable à partir de [RGT] Lot	benu	Target Sollwert Valor asignado Cible	SD SD SD  ET	Range Bereich Intervalo  Intervalle	Target Sollwert Valor asignado Cible	SD SD SD  ET	Range Bereich Intervalo  Intervalle
				Level 1			Level 2		
84600	TSH	20001	µIU/ml	4.35	0.57	2.64 - 6.06	34.95	3.25	25.20 - 44.70
84615	T3	20001	ng/ml	1.16	0.15	0.70 - 1.62	4.48	0.59	2.72 - 6.24
84620	ft3	20001	pg/ml	3.57	0.47	2.16 - 4.98	11.94	1.11	8.61 - 15.27
84610	ft4	22001	pg/ml	10.31	1.35	6.25 - 14.37	34.63	3.22	24.97 - 44.29
82120	FSH	22002	mIU/ml	4.98	0.65	3.02 - 6.94	52.32	4.87	37.72 - 66.92
82115	LH	22002	mIU/ml	4.48	0.59	2.72 - 6.24	61.97	5.76	44.68 - 79.26
84605	T4	20001	ng/ml	63.24	5.88	45.60 - 80.88	156.13	14.52	112.57 - 199.69
84625	anti-TG	22002	IU/ml	32.51	3.02	23.44 - 41.58	102.13	9.50	73.64 - 130.62
84630	anti-TPO	21001	U/ml	16.76	1.56	12.08 - 21.44	28.61	2.66	20.63 - 36.59
83100	AFP	22002	ng/ml	7.13	0.94	4.32 - 9.94	130.40	12.13	94.02 - 166.78
82100	HCG	22002	mIU/ml	2.55	0.33	1.55 - 3.55	79.14	7.36	57.06 - 101.22
85000	25-OH Vitamin D	21001	ng/ml	38.75	3.60	27.94 - 49.56	44.37	4.13	31.99 - 56.75
85055	Vitamin B12	22002	pg/ml	470.82	43.79	339.46 - 602.18	1017.85	94.66	733.87 - 1301.83
82105	Progesterone	21002	ng/ml	0.99	0.13	0.59 - 1.39	22.92	2.13	16.53 - 29.31
83125	fPSA**	21001	ng/ml	1.19	0.16	0.72 - 1.66	3.53	0.46	2.14 - 4.92
83130	PSA**	21001	ng/ml	5.82	0.76	3.53 - 8.11	18.48	1.72	13.32 - 23.64
82125	Prolactin	20001	ng/ml	13.05	1.71	7.91 - 18.19	21.30	1.98	15.36 - 27.24
82110	Testosterone	21001	ng/ml	7.59	1.00	4.60 - 10.58	1.86	0.24	1.13 - 2.59
85030	Troponin I	21001	ng/ml	0.049	0.037	0.00 - 0.159	7.848	0.730	5.658 - 10.038
85040	CK-MB	21001	ng/ml	4.72	0.62	2.86 - 6.58	14.53	1.35	10.48 - 18.58
85035	Myoglobin	21001	ng/ml	39.57	1.96	33.69 - 45.44	176.06	9.49	147.61 - 204.52
83110	CA 125	21001	U/ml	40.43	3.76	29.15 - 51.71	128.73	11.97	92.81 - 164.65
83115	CA 15-3	21001	U/ml	27.24	2.53	19.64 - 34.84	51.10	4.75	36.84 - 65.36
83120	CA 19-9	for 21001 only	U/ml	23.07	2.15	16.63 - 29.51	73.32	6.82	52.86 - 93.78
83120	CA 19-9	22001	U/ml	19.07	1.77	13.75 - 24.39	59.03	5.49	42.56 - 75.50
82130	Estradiol	22002	pg/ml	81.55	8.97	54.64 - 108.46	254.97	28.05	170.83 - 339.11
83105	CEA	22001	ng/ml	2.45	0.32	1.49 - 3.41	9.72	1.28	5.89 - 13.55
85010	PTH	22001	pg/ml	29.05	2.70	20.95 - 37.15	199.76	18.58	144.03 - 255.49

\*\* No CE – only for use outside Europe

