


# HbA1c liquidirect Controls

**REF** 10775

**Content** 2 x 0.5 ml

**LOT** 21002  
 29.02.2024

Version 1 / 26-07-2021

## Target Values / Sollwerte / Valores Asignados/ Valeurs Cibles

Unit
Einheit
Unidad
Unité

Control Blood normal <b>CBN</b>	
Mean	Range
Mittelwert	Bereich
Media	Rango
Moyenne	Plage de valeurs

Control Blood abnormal <b>CBA</b>	
Mean	Range
Mittelwert	Bereich
Media	Rango
Moyenne	Plage de valeurs

Device Gerät Instrumento Equipement	HumaStar 300SR	mmol/mol Hb
		%
	HumaStar 600	mmol/mol Hb
		%
other	mmol/mol Hb	
	%	

38,8	34,0	—	43,5
5,64	5,21	—	6,07
46,5	45,7	—	47,3
6,34	6,26	—	6,41
39,3	32,8	—	45,8
5,69	5,10	—	6,27

118	112	—	125
12,8	12,2	—	13,4
124	121	—	128
13,3	13,0	—	13,7
118	112	—	125
12,8	12,2	—	13,4