

# HUMAN

Target Values

## BioRad Lyphocheck Immunoassay Plus Control

**REF** 370, **LOT** 40350 Level 1 - 3  2020-05-31

REF		Units Einheiten Unitades Unités	Target Sollwert Valor meta Valeur souhaitée	Range Bereich Rango Marge	Target Sollwert Valor meta Valeur souhaitée	Range Bereich Rango Marge	Target Sollwert Valor meta Valeur souhaitée	Range Bereich Rango Marge
<b>HUMAN Hormones &amp; Tumor Markers</b>								
			Level 1 - 40351		Level 2 - 40352		Level 3 - 40353	
52010	AFP	ng/ml	15.6	10.4 - 20.7	89.3	59.8 - 119	203	136 - 270
52020	CEA	ng/ml	1.90	1.27 - 2.53	12.7	8.53 - 16.9	20.4	13.7 - 27.2
55050	Cortisol	ng/ml	41.9	26.8 - 57.0	227	145 - 309	354	226 - 481
55060	DHEA-S	µg/ml	0.59	0.32 - 0.86	1.36	0.75 - 1.98	4.07	2.24 - 5.90
55030	Estradiol	pg/ml	273	153 - 392	713	399 - 1026	1105	619 - 1591
55040	Estriol	ng/ml	1.74	0.95 - 2.52	5.69	3.13 - 8.25	13.6	7.45 - 19.6
52035	fPSA	ng/ml	0.46	0.31 - 0.61	2.16	1.45 - 2.87	>10	-
53020	FSH	IU/l	5.9	3.95 - 7.85	22.0	14.7 - 29.3	46.2	31.0 - 61.5
54015	fT3	pg/ml	2.42	1.62 - 3.22	4.74	3.18 - 6.30	7.08	4.74 - 9.42
54025	fT4	ng/dl	0.95	0.64 - 1.26	2.00	1.34 - 2.66	3.80	2.55 - 5.05
53040	hCG	IU/l	6.71	4.50 - 8.92	20.8	13.9 - 27.7	175	117 - 233
53010	LH	IU/l	1.91	1.28 - 2.54	24.1	16.2 - 32.1	68.0	45.6 - 90.4
53030	Progesterone	ng/ml	1.07	0.59 - 1.55	12.5	6.89 - 18.2	36.8	20.2 - 53.4
55020	Prolactin	ng/ml	8.27	5.54 - 11.0	23.1	15.5 - 30.8	50.7	33.9 - 67.4
52030	PSA	ng/ml	1.20	0.66 - 1.74	4.88	3.17 - 6.59	28.5	18.6 - 38.5
54010	T3	ng/ml	0.80	0.54 - 1.06	2.01	1.35 - 2.67	3.10	2.08 - 4.12
54020	T4	µg/dl	5.50	3.69 - 7.32	9.50	6.37 - 12.6	14.8	9.92 - 19.7
55010	Testosterone	ng/ml	1.92	1.11 - 2.73	6.78	3.93 - 9.62	12.4	7.20 - 17.6
54030	TSH	mIU/l	0.37	0.20 - 0.53	5.26	3.42 - 7.10	*	-

\* Control is not suitable for this parameter as the target value is out of the linearity range of the test.

