

# HUMAN

Target Values

## BioRad Lyphocheck Immunoassay Plus Control

**REF** 370, **LOT** 40380 Level 1 - 3 📅 2022-08-31

Revision 1

REF		Units Einheiten Unitades Unités	Target Sollwert Valor meta Valeur souhaitée	Range Bereich Rango Marge	Target Sollwert Valor meta Valeur souhaitée	Range Bereich Rango Marge	Target Sollwert Valor meta Valeur souhaitée	Range Bereich Rango Marge
<b>HUMAN Hormones &amp; Tumor Markers</b>								
			Level 1 - 40381		Level 2 - 40382		Level 3 - 40383	
52010	AFP	ng/ml	16.9	11.3 - 22.5	100.4	67.2 - 133.5	226.0	151.4 - 300.5
52020	CEA	ng/ml	1.41	0.94 - 1.87	13.1	8.78 - 17.42	21.9	14.7 - 29.1
55050	Cortisol	ng/ml	40.4	25.8 - 54.9	211	135 - 287	345	221 - 469
55060	DHEA-S	µg/ml	0.60	0.33 - 0.86	1.75	0.96 - 2.53	5.43	2.98 - 7.87
55030	Estradiol	pg/ml	170	95.4 - 245	490	274 - 705	922	516 - 1328
55040	Estriol	ng/ml	1.70	0.93 - 2.46	5.12	2.82 - 7.42	12.0	6.63 - 17.5
52035	fPSA	ng/ml	0.43	0.29 - 0.57	2.61	1.75 - 3.47	>10	- - -
53020	FSH	IU/l	5.19	3.48 - 6.91	21.7	14.6 - 28.9	45.7	30.6 - 60.8
54015	ft3	pg/ml	2.35	1.57 - 3.12	5.48	3.67 - 7.29	7.75	5.19 - 10.3
54025	ft4	ng/dl	1.00	0.67 - 1.33	2.35	1.58 - 3.13	3.73	2.50 - 4.96
53040	hCG	IU/l	6.37	4.27 - 8.48	24.6	16.5 - 32.7	173.3	116.1 - 230.4
53010	LH	IU/l	1.89	1.26 - 2.51	29.6	19.8 - 39.4	93.1	62.4 - 123.8
53030	Progesterone	ng/ml	0.95	0.52 - 1.37	14.5	7.99 - 21.1	<40	- - -
55020	Prolactin	ng/ml	7.60	5.09 - 10.1	22.0	14.7 - 29.2	45.5	30.5 - 60.5
52030	PSA	ng/ml	1.49	0.82 - 2.16	5.11	3.32 - 6.90	29.8	19.4 - 40.2
54010	T3	ng/ml	0.91	0.55 - 1.27	2.47	1.65 - 3.28	3.73	2.50 - 4.96
54020	T4	µg/dl	5.54	3.71 - 7.36	10.2	6.84 - 13.6	14	9.38 - 18.6
55010	Testosterone	ng/ml	1.57	0.91 - 2.23	4.65	2.70 - 6.61	10.5	6.12 - 15.0
54030	TSH	mIU/l	0.29	0.16 - 0.41	5.31	3.45 - 7.17	28.8	18.7 - >30

\* The marked level is not applicable for the indicated parameter as the target value is out of the measuring range of the assay.

