

HUMAN

Target Values

BioRad Lyphocheck Immunoassay Plus Control

370 Trilevel 12x5ml
371 Level 1 12x5ml
372 Level 2 12x5ml
373 Level 3 12x5ml
370X MiniPak 3x5ml

Version 04 / 05-2022

REF

LOT

40390

Level 1 40391

Level 2 40392

Level 3 40393

IVD



2022-12-31

Link for latest version of target value sheet:

<https://www.human.de/clia/target-value-sheets/>



REF	Assay Parameter Ensayo Essai	Valid from Gültig ab Válido a partir de Valable à partir de RGT Lot	Units Einheiten Unidades Unité	Target Sollwert Valor asignado Cible	Range Bereich Intervalo Intervalle	Target Sollwert Valor asignado Cible	Range Bereich Intervalo Intervalle	Target Sollwert Valor asignado Cible	Range Bereich Intervalo Intervalle
				Level 1 - 40391		Level 2 - 40392		Level 3 - 40393	
84600	TSH	20001	µIU/ml	0.660	0.430 - 0.900	5.89	3.87 - 7.99	37.3	24.5 - 50.6
84620	ft3	20001	pg/ml	2.50	1.64 - 3.39	5.93	3.90 - 8.05	8.46	5.56 - 11.5
84610	ft4	20001	pg/ml	5.83	3.83 - 7.91	13.0	8.55 - 17.6	19.8	13.0 - 26.9
82120	FSH	20001	mIU/ml	5.50	3.62 - 7.46	22.8	15.0 - 30.9	51.8	34.1 - 70.3
84605	T4	20001	ng/ml	50.0	32.9 - 67.8	94.3	62.0 - 128	112	73.3 - 151
83130	PSA	20001	ng/ml	1.06	0.700 - 1.44	4.54	2.98 - 6.16	27.0	17.7 - 36.6
84615	T3	20001	ng/ml	0.81	0.55 - 1.07	1.98	1.45 - 2.51	2.99	2.18 - 3.80
82115	LH	20001	mIU/ml	1.88	1.31 - 2.44	26.43	18.39 - 34.31	103.33	71.88 - 134.14
83125	fPSA	20001	ng/ml	0.29	0.19 - 0.39	1.85	1.30 - 2.41	10.97	7.68 - 14.26
83100	AFP	20001	ng/ml	18.26	13.33 - 23.27	105.61	77.10 - 134.59	227.25	165.89 - 289.60
82100	HCG	20001	mIU/ml	4.91	3.27 - 6.54	25.12	16.73 - 33.45	163.69	108.99 - 217.97
82105	Progesterone	20001	ng/ml	0.27	0.18 - 0.36	11.91	8.34 - 15.45	31.21	21.86 - 40.49
82110	Testosterone	20001	ng/ml	0.71	0.48 - 0.95	3.34	2.22 - 4.45	7.71	5.14 - 10.27
82125	Prolactin	20001	ng/ml	9.22	6.46 - 11.98	25.96	19.42 - 32.50	54.86	41.04 - 68.68
84625	anti-TG	20001	IU/ml	101.66	75.93 - 126.76	94.82	70.82 - 118.23	67.94	50.75 - 84.72
84630	anti-TPO	20001	U/ml	29.40	21.96 - 36.66	24.93	18.62 - 31.09	21.78	16.27 - 27.16
85000	25-OH Vitamin D	20001	ng/ml	21.05	15.72 - 26.25	28.60	21.36 - 35.66	23.79	17.77 - 29.66
83105	CEA	22001	ng/ml	2.08	1.37 - 2.82	14.5	9.5 - 19.6	38.6	25.4 - 52.4
82130	Estradiol	22001	pg/ml	53.5	34.8 - 72.2	198	128 - 267	575	374 - 775