

REF 13010

Target Values / Sollwerte / Valores Meta/ Valeurs Souhaitées

Please note: temporary target values of RF are valid until 30.11.2019

The latest version of the RF target value sheet is accessible via:

<https://www.human.de/target-value-sheets/>

HUMAN Constituent / Method		SI Units SI Einheiten Unitades SI Unités SI	Target Sollwert Valor meta Valeur souhaitée	Range Bereich Rango Marge	Units Einheiten Unitades Unités	Target Sollwert Valor meta Valeur souhaitée	Range Bereich Rango Marge
Level 1							
RHEUMATOID FACTORS (RF) Immunoturbidimetry	Multi purpose reagent HS 600 (system reagent)	klU/l klU/l	25.4 25.4	20.3 - 30.5 20.3 - 30.5	IU/ml IU/ml	25.4 25.4	20.3 - 30.5 20.3 - 30.5
RHEUMATOID FACTORS (RF) Immunoturbidimetry latex enhanced	Multi purpose reagent	klU/l	26.8	21.4 - 32.2	IU/ml	26.8	21.4 - 32.2
Level 2							
RHEUMATOID FACTORS (RF) Immunoturbidimetry	Multi purpose reagent HS 600 (system reagent)	klU/l klU/l	81.1 81.1	64.9 - 97.3 64.9 - 97.3	IU/ml IU/ml	81.1 81.1	64.9 - 97.3 64.9 - 97.3
RHEUMATOID FACTORS (RF) Immunoturbidimetry latex enhanced	Multi purpose reagent	klU/l	82.8	66.2 - 99.4	IU/ml	82.8	66.2 - 99.4