

TURBIDOS

Target Values / Sollwerte / Valores Meta/ Valeurs Souhaitées

IVD
Lot 17002

REF 13010

Please note: temporary target values of RF are valid until 30.04.2018.

The latest version of the RF target value sheet is accessible via:

www.human.de/data/gb/vr/targetvaluerf-17002.pdf

www.human-de.com/data/gb/vr/targetvaluerf-17002.pdf

HUMAN Constituent / Method		SI Units SI Einheiten Unitades SI Unités SI	Target Sollwert Valor meta Valeur souhaitée	Range Bereich Rango Marge	Units Einheiten Unitades Unités	Target Sollwert Valor meta Valeur souhaitée	Range Bereich Rango Marge
Level 1							
RHEUMATOID FACTORS (RF)	Multi purpose reagent	klU/l	18,7	14,9 - 22,4	IU/ml	18,7	14,9 - 22,4
Immunoturbidimetry	HumaStar 600 system reagent	klU/l	20,7	16,6 - 24,9	IU/ml	20,7	16,6 - 24,9
RHEUMATOID FACTORS (RF)	Multi purpose reagent	klU/l	27,1	21,7 - 32,6	IU/ml	27,1	21,7 - 32,6
Immunoturbidimetry latex enhanced							
Level 2							
RHEUMATOID FACTORS (RF)	Multi purpose reagent	klU/l	71,1	56,8 - 85,3	IU/ml	71,1	56,8 - 85,3
Immunoturbidimetry	HumaStar 600 system reagent	klU/l	61,4	49,1 - 73,7	IU/ml	61,4	49,1 - 73,7
RHEUMATOID FACTORS (RF)	Multi purpose reagent	klU/l	86,0	68,8 - 103,2	IU/ml	86,0	68,8 - 103,2
Immunoturbidimetry latex enhanced							

CS-TU-17002 INF 1301004 12-2017-03



Human Gesellschaft für Biochemica und Diagnostica mbH
Max-Planck-Ring 21 · 65205 Wiesbaden · Germany