

# HUMAN

Target Values

## BioRad Lyphocheck Tumor Marker Plus Control

Version 02/ 05 - 2022

367 Level 1 6x2ml

Level 1 54691


**REF**

368 Level 2 6x2ml

**LOT** 54690

Level 2 54692

**IVD**

 2023-06-30

369 Level 3 6x2ml

Level 3 54693

368X Trilevel MiniPak 3x2ml

Link for latest version of target value sheet:

<https://www.human.de/clia/target-value-sheets/>



<b>REF</b>	Assay Parameter Ensayo  Essai	Valid from Gültig ab Válido a partir de  Valable à partir de [RGT] Lot	Units Einheiten Unidades  Unité	Target Sollwert Valor asignado Cible	Range Bereich Intervalo  Intervalle	Target Sollwert Valor asignado Cible	Range Bereich Intervalo  Intervalle	Target Sollwert Valor asignado Cible	Range Bereich Intervalo  Intervalle
				Level 1 - 54691		Level 2 - 54692		Level 3 - 54693	
83130	PSA	20001	ng/ml	0.07	0.01 - 0.13	3.22	2.12 - 4.37	15.78	10.37 - 21.42
83125	fPSA	20001	ng/ml	0.03	< 0,01 - 0.06	2.75	1.81 - 3.73	12.66	8.32 - 17.18
83100	AFP	20001	ng/ml	9.66	6.47 - 12.85	74.04	49.61 - 98.47	207.49	139.02 - 275.96
83115	CA 15-3	21001	U/ml	18.5	12.2 - 25.1	46.4	30.5 - 63.0	97.9	64.3 - 133.0
83105	CEA	22001	ng/ml	2.65	1.74 - 3.60	18.9	12.4 - 25.7	60.6	39.8 - 82.2

#54690

Rev. 02 | valid of 13.05.2022



Human Gesellschaft für Biochemica und Diagnostica mbH

Max-Planck-Ring 21 - D-65205 Wiesbaden - Germany

Tel.: +496122 99880 Fax: +496122 9988100 - e-mail: human@human.de